



# Technical Information

## Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes



### OVERVIEW

Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes are engineered as ready-to-use cleaner-disinfectant wipes for healthcare facilities. They are compatible with a wide range of common healthcare surfaces and have a low residue and odor profile.

These EPA-registered disinfectant wipes contain hydrogen peroxide and other ingredients to kill most bacteria and viruses on a variety of hard, nonporous healthcare surfaces with contact times of 30 seconds or 1 minute.

### INGREDIENT INFORMATION

Active ingredient: Hydrogen Peroxide..... 1.4%  
 Other ingredients:..... 98.6%  
 Total (not including weight of wipe):..... 100.00%

### AVAILABLE FORMATS

Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes are available in multiple sizes and formats to accomplish different disinfecting jobs.

Format	Image	Designed for	Item #	Quantity	Shelf Life
<b>Clinical Surface Wipes</b> (6.75" x 5.75")		Daily Disinfection Small high-touch surfaces Medical equipment	30825	155 ct. canister, 6/case	24 months
<b>Multi-Purpose Wipes</b> (6.75" x 9")		Daily Disinfection Patient room surfaces Medical equipment	30824	95 ct. canister, 6/case	24 months

At a Glance Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes	
EPA Registration Number	7619-25
Product Classification	Pre-moistened wipes
Dilution	Ready to use
Disinfectant Type	1-step cleaner-disinfectant
Active Ingredient Chemistry	Hydrogen Peroxide
Odor	Cleaning Agent
Shelf Life	Expiration date listed on each package: month, day, year
Toxicity Signal Word	CAUTION
HMIS Rating	1-0-0-A
PPE Required	Follow your facility's protocols



# Technical Information

## Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes



PRODUCTS  
& SAFETY

### DISINFECTION EFFICACY CLAIMS FOR WIPES

	Disinfecting efficacy against:	KILL TIMES
Bacteria	<i>Acinetobacter baumannii</i>	1 min
	<i>Campylobacter jejuni</i>	1 min
	<i>Enterobacter aerogenes</i>	30 sec
	Vancomycin-resistant <i>Enterococcus faecalis</i> (VRE)	30 sec
	Multidrug-resistant <i>Enterococcus faecium</i>	30 sec
	ESBL-producing <i>Escherichia coli</i>	30 sec
	<i>Escherichia coli</i>	30 sec
	<i>Escherichia coli</i> O157:H7	30 sec
	<i>Escherichia coli</i> (NDM-1)	30 sec
	<i>Klebsiella oxytoca</i>	30 sec
	Carbapenem-resistant <i>Klebsiella pneumoniae</i>	30 sec
	<i>Klebsiella pneumoniae</i> (NDM-1)	30 sec
	ESBL-producing <i>Klebsiella pneumoniae</i>	30 sec
	Multidrug-resistant <i>Klebsiella pneumoniae</i>	30 sec
	<i>Pseudomonas aeruginosa</i>	30 sec
	<i>Salmonella enterica</i>	30 sec
	<i>Staphylococcus aureus</i>	1 min
	Linezolid-resistant <i>Staphylococcus aureus</i> (LRSA)	1 min
	Methicillin-resistant <i>Staphylococcus aureus</i>	1 min
	Methicillin-resistant <i>Staphylococcus aureus</i> (CA-MRSA USA300)	1 min
<i>Stenotrophomonas maltophilia</i>	30 sec	
Penicillin-resistant <i>Streptococcus pneumoniae</i>	30 sec	
<i>Streptococcus pyogenes</i>	30 sec	
Viruses	Avian Influenza A Virus (H5N1)	1 min
	Cytomegalovirus	30 sec
	Hepatitis B Virus (HBV) <sup>†</sup>	30 sec
	Hepatitis C Virus (HCV)	30 sec
	Herpes Simplex Virus Type 1 (HSV-1)	30 sec
	Herpes Simplex Virus Type 2 (HSV-2)	30 sec
	Human Immunodeficiency Virus Type 1 (HIV-1)	30 sec
	Human Coronavirus	30 sec
	Influenza A Virus (H1N1)	30 sec
	Influenza A Virus/Hong Kong	30 sec
	Influenza B Virus	30 sec
	Norovirus <sup>‡</sup>	3 min
	Respiratory Syncytial Virus (RSV)	30 sec
	Rhinovirus	1 min
	Rotavirus	1 min
Fungi	<i>Candida albicans</i>	5 min
	<i>Trichophyton mentagrophytes</i> (Athlete's Foot Fungus)	5 min
TB	<i>Mycobacterium bovis</i> (TB)	5 min

<sup>†</sup> Obtained via a surrogate using feline calicivirus.

<sup>‡</sup> [www.osha.gov/SLTC/healthguidelines/isopropylalcohol/recognition.html](http://www.osha.gov/SLTC/healthguidelines/isopropylalcohol/recognition.html)

Use as directed on hard, nonporous surfaces. All claims based on data as of 2/2014.



## Technical Information

### Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes



PRODUCTS  
& SAFETY

## SAFETY INFORMATION

Reference the Safety Data Sheet (SDS) sheet for complete product safety information.

### EPA Label Statements

This chemical is a pesticide product registered by the Environmental Protection Agency and is subject to certain labeling requirements under federal pesticide law. These requirements differ from the classification criteria and hazard information required for safety data sheets and for workplace labels of non-pesticide chemicals. Following is the hazard information as required on the pesticide label:

**CAUTION:** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

**FIRST AID: IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or when going for treatment.

**PHYSICAL AND CHEMICAL HAZARDS:** The product contains no substances which at their given concentration, are considered to be hazardous to health

### National Fire Protection Association (NFPA) Ratings

Health Hazard: 1

Flammability: 0

Instability: 0

Physical and Chemical Hazards: None

### Hazardous Materials Identification System (HMIS) Ratings

Health Hazard: 1

Flammability: 0

Physical Hazard: 0

Personal Protection: A

Questions?  
Call (800) 243-5350

Distributed By:



Contact your Pharmacal sales representative  
at 800-243-5350 for further information.

Use products as directed

© 2017 Clorox Professional Products Company