

# CLIDOX-S<sup>®</sup> MIXING STATION

Effectively measures and mixes CLIDOX-S<sup>®</sup> automatically.



Precise Measurements

Reliable Quality

Laboratory Results

CLIDOX-S<sup>®</sup> Mixing Station allows for total product quality assurance and safety. CLIDOX-S<sup>®</sup> (Chlorine Dioxide) must be prepared on-site, at the correct dilution rates, and in the proper sequence. CLIDOX-S<sup>®</sup> Mixing Station performs these functions flawlessly and accurately.

CLIDOX-S<sup>®</sup> Mixing Station can be custom designed to:

- prepare CLIDOX-S<sup>®</sup> at 1:3:1, 1:5:1, or 1:18:1 use dilution
- fill various size container
- have multiple programs for varying needs
- have low level alarms for the chemical supply
- have an automatic fresh water fill

For the critical difference, call Pharmacal for an evaluation of your needs, options and price quotation.

**pharmacal**  
RESEARCH LABS. INC.

562 Captain Neville Dr. Waterbury CT. 06705  
[www.pharmacal.com](http://www.pharmacal.com) 800-243-5350